

# Water Analysis

VISOCOLOR® SCHOOL



Water analysis - it's that simple

## VISOCOLOR® SCHOOL reagent case

- Compact mini lab
- Easy test procedures
- Safe results
- Flexible applications
- Incredible value

MACHEREY-NAGEL

[www.mn-net.com](http://www.mn-net.com)

Distribútor pre SR:  
AZ CHROM s.r.o., Robotnícka 10, 831 03 Bratislava  
Tel. 0907 244526 [azetchrom@hplc.sk](mailto:azetchrom@hplc.sk), [www.azetchrom.sk](http://www.azetchrom.sk)



# VISOCOLOR® SCHOOL

## Water analysis for kids

### Compact reagent case

The **VISOCOLOR® SCHOOL** reagent case is especially designed for schools and caters to the needs of both students and teachers. All reagents are approved to be used in schools in Germany (GUV-SR 2004 directive) and can be disposed of easily just down the drain without any harm to the environment. The case contains 6 colorimetric and titrimetric tests for at least 50 determinations each to determine the most important water parameters:

<b>Ammonium</b>	0.2–3 mg/L $\text{NH}_4^+$
<b>Nitrate</b>	1–90 mg/L $\text{NO}_3^-$
<b>Nitrite</b>	0.02–0.5 mg/L $\text{NO}_2^-$
<b>Phosphate</b>	0.5–15 mg/L $\text{PO}_4^{3-}$
<b>pH</b>	4.0–9.0
<b>Total hardness</b>	1 drop = 17.8 mg/L $\text{CaCO}_3$




The newly developed tests are based on our proven **VISOCOLOR®** system and were developed as high quality **VISOCOLOR® SCHOOL** test specifically for schools.


## Your benefits

### Great features – incredible price

- Maximum safety due to exact labeling of all reagents
- Safe results using color and turbidity compensation
- Especially stable and rugged case as well as chemical resistant foam inlay
- High sensitivity down to the values of German drinking water standards
- Safe for the environment and easy disposal of used tests



Precise color matching with all-in-one color scale



Easy test procedures with pictogram instructions



Extensive manual with valuable background information

Optionally, the case can be equipped with a pack of pH-Fix or QUANTOFIX® test strips as well.

Thus, the **VISOCOLOR® SCHOOL** reagent case is the ideal tool for projects revolving around water and water analysis in any science class.

### Ordering information

#### **VISOCOLOR® SCHOOL** reagent case

including all necessary reagents and accessories for ammonium, total hardness, nitrate, nitrite, phosphate and pH, manual and color scale

REF 933 100

Distributed by:

#### **VISOCOLOR® SCHOOL** refill pack

REF 933 200

#### **VISOCOLOR® SCHOOL** color chart

REF 933 300

#### **VISOCOLOR® SCHOOL** manual (German / English)

REF 933 150

# MACHEREY-NAGEL



Distribútor pre SR:  
AZ CHROM s.r.o., Robotnícka 10, 831 03 Bratislava  
Tel. 0907 244526 azetchrom@hplc.sk, www.azetchrom.sk

